



CERTIFICATE OF ANALYSIS

Solution Tampon pH 4.01 ± 0.02 (25 °C)
Selon NIST

Reference	CLPH401.1K	Production date	mai-25
Batch N°	CMPH401.V0525	Expiry date	juin-27
Grade	Buffer Solution	Molecular formula	-
CAS Number	-	Molecular mass	-

Specifications	Unit	Standards	Results
pH (25°C) (tolérance ± 0.02)	pH	3.99 - 4.03	4.01
pH laboratory uncertainty	pH	± 0.02 (k=2; 95 %)	± 0.02 (k=2; 95 %)

pH	T(°C)	pH	T(°C)
4.01	0	4.00	20
4.00	5	4.01	25
4.00	10	4.02	30
4.00	15	4.03	35

Additional Information:

This product is conform to the above specifications

The document has been produced electronically and is valid without signature.

CLCM LABS
Quality Management

