



CERTIFICATE OF ANALYSIS

Acetone
≥ 98% Extra Pure
Accord. Ph. Eur.

Reference	CL5808.2K5	Production date	Jun-25
Batch N°	CM5808.D0625	Expiry date	Jun-28
Grade	Lab reagent	Molecular formula	C ₃ H ₆ O
CAS Number	67-64-1	Molecular mass	58,08 g/mol

Specifications	Unit	Standards	Results
Purity	%	99,5 min	99.7
Appearance	Visual	Clean and clear	Complies
Color	Hazen	10 max	5
Water	%	0,3 max	0.22
Refractive index at 20 °C	/	Not specified(ASTM D 1218)	Complies
Density at 20 °C	g/ml	0,789 - 0,798	0.791
Acidity	ppm	20 max	< 20

Additional Information :

This product is conform to the above specifications

The document has been produced electronically and is valid without signature.

CLCM LABS

Quality Management

