



CERTIFICATE OF ANALYSIS

Formaldehyde
Méthanal
Stabilised with methanol 36%

Reference	CL3003.1K	Production date	Feb-25
Batch N°	CM3003.D022536	Expiry date	May-28
Grade	Lab Reagent	Molecular formula	CH ₂ O
CAS Number	50-00-0	Molecular mass	30.03 g/mol

Specifications	Unit	Standards	Results
Assay	%	36.0-38.0	Complies
Identity	-	Must comply	Complies
Appearance of solution	-	Must comply	Complies
Sulphated ash	%	≤ 0.005	< 0.005
Acidity	-	Must comply	Complies
Methanol (GC)	%	9.0-15.0	10
Chloride (Cl)	%	≤ 0.003	< 0.003
Sulphate (SO ₄)	%	≤ 0.005	< 0.005

Additional Information :

This product is conform to the above specifications

The document has been produced electronically and is valid without signature.

CLCM LABS
Quality Management

