



Product of european origin conditioned by CLCM-LABS
SARL CLCM-labs City Sadikia hangar N°04,
Laghouat 03000 - Algérie
Fax : 029113704 / num : 0555655655
contact@clcm lab.dz



CERTIFICATE OF ANALYSIS

Conductivity Standard solution

Solution standard de conductivité

12880 $\mu\text{S}/\text{cm}$ (25 °C)

Selon NIST

Reference
Batch N°
Grade

CLEC12880.500
CMEC12880.V0525
Lab Standard

Production date
Expiry date

May-25
May-27

Specifications	Unit	Standards	Results
Conductivity @25°C	$\mu\text{S}/\text{cm}$	12880 $\pm$ 1%	12643
Conductivity @18°C	$\mu\text{S}/\text{cm}$	11174 $\pm$ 1%	10890
Appearance	Visual	Clear colorless liquide	Complies
Odor	/	Odorless	Complies
Density @25°C	g/ml	1.00 - 1.01	1.00
Refractive index	Brix	1.332 - 1.334	1.334
Conductivity.($\mu\text{S}/\text{cm}$)	T(°C)	Conductivity.($\mu\text{S}/\text{cm}$)	T(°C)
10439	15	11909	21
10648	16	12153	22
10929	17	12358	23
11174	18	12643	24
11419	19	12880	25
11664	20	14112	30

Additional Information :

This product is conform to the above specifications

The document has been produced electronically and is valid without signature.

CLCM LABS
Quality Management



SARL CLCM-labs, Cité Sadikia hangar N°04, Laghouat 03000 - Algérie