



Product of european origin conditioned by CLCM-LABS
SARL CLCM-labs City Sadikia hangar N°04,
Laghouat 03000 - Algérie
Fax : 029113704 / num : 0555655655
contact@clcm lab.dz



CERTIFICATE OF ANALYSIS

Conductivity Standard solution

Solution standard de conductivité

1413 $\mu\text{S}/\text{cm}$ (25 °C)

Selon NIST

Reference
Batch N°
Grade

CLEC1413.500
CMEC1413.V0525
Lab Standard

Production date
Expiry date

May-25
May-27

Specifications	Unit	Standards	Results
Conductivity @25°C	$\mu\text{S}/\text{cm}$	1413 $\pm$ 1%	1407
Conductivity @18°C	$\mu\text{S}/\text{cm}$	1221 $\pm$ 1%	1218
Appearance	Visual	Clear colorless liquide	Complies
Odor	/	Odorless	Complies
Density @25°C	g/ml	1.00 - 1.01	1.00
Refractive index	Brix	1.332 - 1.334	1.334
Conductivity.($\mu\text{S}/\text{cm}$)	T(°C)	Conductivity.($\mu\text{S}/\text{cm}$)	T(°C)
1139	15	1304	21
1167	16	1331	22
1194	17	1358	23
1221	18	1386	24
1249	19	1413	25
1276	20	1550	30

Additional Information :

This product is conform to the above specifications

The document has been produced electronically and is valid without signature.

CLCM LABS

Quality Management



SARL CLCM-labs, Cité Sadikia hangar N°04, Laghouat 03000 - Algérie