



Product of european origin conditioned by CLCM-LABS
SARL CLCM-labs City Sadikia hangar N°04,
Laghouat 03000 - Algérie
Fax : 029113704 / num : 0555655655
contact@clcm lab.dz



CERTIFICATE OF ANALYSIS

Conductivity Standard solution

Solution standard de conductivité

147 $\mu\text{S}/\text{cm}$ (25 °C)

Selon NIST

Reference
Batch N°
Grade

CLEC147.500
CMEC147.V0925
Lab Standard

Production date
Expiry date

sept-25
déc-27

Specifications	Unit	Standards	Results
Conductivity @25°C	$\mu\text{S}/\text{cm}$	147 $\pm$ 1%	146,75
Conductivity @18°C	$\mu\text{S}/\text{cm}$	126.91 $\pm$ 1%	126.70
Appearance	Visual	Clear colorless liquide	Complies
Odor	/	Odorless	Complies
Density @25°C	g/ml	1.00 - 1.01	1.00
Refractive index	Brix	1.332 - 1.334	1,334
Conductivity.($\mu\text{S}/\text{cm}$)	T(°C)	Conductivity.($\mu\text{S}/\text{cm}$)	T(°C)
118,31	15	135,52	21
121,18	16	138,39	22
124,05	17	141,26	23
126,91	18	144,13	24
129,78	19	147	25
132,69	20	161,33	30

Additional Information :

This product is conform to the above specifications

The document has been produced electronically and is valid without signature.

CLCM LABS
Quality Management



SARL CLCM-labs, Cité Sadikia hangar N°04, Laghouat 03000 - Algérie