

SLY-C PLUS AUTOMATIC DIGITAL SEED COUNTER

Automatic seed Counter can automatic count many kinds of grains by setting number. In agriculture it is mainly used in grain count, exp. Rice, wheat, sorghum, corn, vegetables seed and so on.

Feature:

- Two modes: Free count and bulk count;
- The counting speed is adjustable;
- High accuracy: Sensitivity levels for different size of particles;
- Automatic deceleration: when counting number is close to the setting number, the counting speed slow down automatically;
- Automatic motor Protection: the motor will stop work once it detects there is no particles on the plate for 60seconds;
- User-friendly interface: LCD screen, mode selection, counting time, level and number of particles display clearly.
- Widely application: suitable for particles includes but is not limited to: crops seeds, vegetable seeds, feed, diamond, tablet.

Technical parameter:

Count Range	1-12mm
Count Accuracy	$\pm 2\%$ $\geq 1000 / 3\text{min}$ (Sesame)
Count Speed	1~99999
Counting Range	1-99999
Bulk Counting	<70db
Noise	AC110V /60Hz and AC 220V /50Hz
Power	$\geq 5\text{Hr}$
Continuous working time	285*210*176mm
Size of Machine	-10°C-50°C, <85%RH
Operation Environment	

